

RESPIRATOR LOCKER®

2250—WITH 2200N95 MASK
2350—WITH 2300N95 MASK

Holds 5 masks per Locker, a week's supply in most work situations*.

Allows employer to meter respirator usage, reducing costs.

Reduces worker time spent going to get new masks.

Handy supply for workers at field worksites.

Provides each worker with hygienic storage and a personal supply of respirators.

Separate chamber stores mask during breaks.



HANDY TWO IN ONE STORAGE!

Respirator Locker provides hygienic storage for masks. It's a personal dispenser for workers required to travel to remote locations. Fits easily in employee lockers or work vehicles.

Separate compartment permits storage of a used mask during breaks, prevents premature disposal of masks.

LESS WASTE, SAVES MONEY.

Reduces employee waste from taking extra masks for the day. Allows employer to meter respirator usage. Saves time going to tool storage for new respirators. The result is fewer respirators used and lower costs.

EASY FIT, MORE COMFORT.

Each Moldex respirator has our exclusive Dura-Mesh® shell that protects the filter

media, so it stays cleaner looking longer. The shell resists collapsing, even in heat and humidity. So the respirator holds its shape longer. Model #2350 also features an exhalation valve that reduces hot air buildup. Workers stay cooler and comfortable longer.

N95 PROTECTION.

NIOSH certified under 42CFR84 regulations to have a filter efficiency level of 95% against particulates free of oil up to 10X PEL.

MOLDEX TECHNICAL SUPPORT.

Moldex provides technical assistance including help with fit testing, a *Train-the-Trainer Program*, a manual on *How to Set Up Your Respiratory Protection Program* plus a Toll-Free Technical Help line. Just call 800/421-0668 ext. 550 or 310/837-6500 ext. 550.

*Masks should be replaced in accordance with all applicable regulations.

**N95
PARTICULATE
RESPIRATORS**

RESPIRATOR LOCKER®

2250—WITH 2200N95 MASK
2350—WITH 2300N95 MASK

APPLICATIONS:

⚠ WARNING TO USER:

Follow all instructions and warnings on the use of this respirator. Failure to do so will reduce respirator effectiveness, wearer protection and may result in sickness or death. Wear during all times of exposure.

Before use, the user must first be trained by the employer in proper respirator use in accordance with applicable safety and health standards for the contaminant and exposure level in the assigned work area.

The particulates and other contaminants which can be dangerous to your health include those which you cannot see, taste, or smell.

This product has not been sold with warnings or use instructions for personnel involved in healthcare or related situations, where there may be the possibility of contact with disease or biological hazards. If you are considering such uses, first call the Moldex Technical Service Dept., 800/421-0668 ext. 550 or 310/837-6500 ext. 550.

RESTRICTIONS:

1. This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
2. Do not remain in contaminated area if any physical distress occurs, for example breathing difficulty, dizziness or nausea.
3. Leave contaminated area and change respirator if it is damaged, distorted, a proper fit can not be obtained or breathing becomes difficult.
4. Do not alter, modify, or abuse this respirator.
5. Dispose of respirator according to your employer's policy and local regulations.
6. Dispose no later than thirty days after first use.
7. Use respirator before the "use by" expiration date printed on box.

USE AGAINST:

1. Particulate aerosols free of oil.

DO NOT USE AGAINST:

1. Concentrations of contaminants which are unknown or are immediately dangerous to life or health.
2. Concentrations of particulates which exceed the maximum use concentration or 10 times the OSHA Permissible Exposure Limit, whichever is lower.
3. Gases, vapors, asbestos, sandblasting, paint spray, or particulate material which generate harmful vapors.
4. Oil-based mists.

NIOSH LABEL



MOLDEX/METRIC INCORPORATED
CULVER CITY, CALIFORNIA, U.S.A.
800/421-0668



THESE RESPIRATORS ARE APPROVED ONLY IN THE FOLLOWING CONFIGURATIONS

TC-	PROTECTION ¹	RESPIRATOR		CAUTIONS AND LIMITATIONS ²
		2200N95	2300N95	
84A-0327	N95	X		A B C J M N O P
84A-0328	N95		X	A B C J M N O P

1. Protection

N95 Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply

2. Cautions and Limitations

A- Not for use in atmospheres containing less than 19.5 percent oxygen.
 B- Not for use in atmospheres immediately dangerous to life or health.
 C- Do not exceed maximum use concentrations established by regulatory standards.
 J- Failure to properly use and maintain this product could result in injury or death.
 M- All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
 N- Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configurations as specified by the manufacturer.
 O- Refer to user's instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
 P- NIOSH does not evaluate respirators as surgical masks.

TECHNICAL HELP LINE

For information, technical assistance, and training materials call 800/421-0668 ext. 550 or 310/837-6500 ext. 550.

PACKAGING

2250	2350
5 Masks per Locker	5 Masks per Locker
27 Lockers per Case	16 Lockers per Case
135 Masks per Case	80 Masks per Case

LIMITED WARRANTY

IMPORTANT NOTICE TO PURCHASER

This limited warranty is made in lieu of the warranties of merchantability, fitness for particular purposes and all other warranties, express or implied. There are no other warranties which extend beyond the description on the face hereof. The physical standards and specifications of Moldex will be met by products sold. **Exclusive Remedies:** damages for the breach of this limited warranty are limited to the replacement of such quantity of Moldex, products proved to be defectively manufactured. Except as provided above, Moldex shall not be liable or responsible for any loss, damage, or liability, direct, indirect, incidental, special, or consequential, arising out of sale, use, or misuse, or the inability to use products by the user.

DISTRIBUTED BY:

In Canada:
Pacesetter Sales & Associates
20 Arthur Hall Drive
Sharon, ON LOG 1V0
Tel: 905/478-8042 Fax: 905/478-8010
www.pacesettersales.com

MOLDEX-METRIC, INC.
10111 W. Jefferson Blvd.
Culver City, CA 90232
800/421-0668 310/837-6500
FAX 310/837-9563

www.moldex.com e-mail: sales@moldex.com
Moldex, Dura-Mesh and Respirator Locker are registered trademarks. The Mesh Design is a trademark of Moldex.
U.S. Patents #D285,374, #D418,410, #4,384,577, #4,850,347, #6,076,698.
Made in the U.S.A. 9701-718 REV D 9/05

