

EZ-ON[®] Series Disposable Respirators



EASY ON, EASY OFF

The first and only NIOSH N95 approved one adjustable cotton strap respirator gives you a quick fit. The single strap with Flexwing[®] and head harness provides an easier and more comfortable fit than two straps. The adjustable cotton strap is cool and comfortable against the skin.

DURA-MESH[®] SHELL LASTS LONGER

Exclusive Dura-Mesh[®] shells resist collapsing, even in heat and humidity so the respirator holds its shape longer. The result is fewer respirators used and lower costs.

VENT HOT AIR, STAY COOLER

Exclusive Ventex[®] valve lets the hot air out faster so breathing resistance is reduced and workers stay cooler.

N95 NON-OIL PROTECTION

NIOSH certified to have a filter efficiency level of 95% or greater against non-oil based particulate aerosols.

APPLICATIONS

Hot, humid and dusty environments under 10 X PEL, i.e. welding, soldering and brazing, grinding, sanding, sweeping, bagging, stone quarrying, metal machining, cement, underground mining, polishing, buffing, lime, poultry and textiles. No oils present.

FEATURES

- 100% PVC-Free.
- The first and only NIOSH N95 approved one adjustable cotton strap respirator.
- Easy on, easy off one strap design for a quick fit.
- Single strap with Flexwing[®] and head harness provides an easier and more comfortable fit than two straps.
- Cotton strap is adjustable for a cool, custom fit.
- Outer Dura-Mesh[®] shell keeps filter media cleaner looking longer.
- Inner Dura-Mesh[®] shell resists collapsing in heat and humidity.
- Exclusive Ventex[®] valve lets the hot air out faster so workers breathe easier and stay cooler.
- Soft foam nose cushion for added comfort.



USE AGAINST:

1. Particulate aerosols free of oil.

DO NOT USE AGAINST:

1. Concentrations of contaminants which are unknown, or are immediately dangerous to life or health.
2. Concentrations of particulates which exceed the maximum use concentration or 10 times the OSHA Permissible Exposure Limit, whichever is lower.
3. Gases, vapors, asbestos, paint spray, sandblasting or particulate materials which generate harmful vapors.
4. Oil-based mists.

RESTRICTIONS

1. This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
2. Do not remain in contaminated area if any physical distress occurs, for example breathing difficulty, dizziness or nausea.
3. Leave contaminated area and replace respirator if it is damaged, distorted, a proper fit can not be obtained or breathing becomes difficult.
4. Prior to each use, carefully inspect the entire respirator, including filter media and strap/weld attachment area for tears and damage. Staple or strap/weld attachment perforations do not affect respirator performance.
5. Do not alter, modify, or abuse this respirator.
6. Store unused respirators in box/bag in a clean, dry, non-contaminated area.
7. Dispose of respirator according to your employer's policy and local regulations.
8. Dispose no later than thirty days after first use.
9. Use respirator before the "use by" expiration date printed on box/bag.
10. If used for welding, wear appropriate eye and face protection.

TECHNICAL SERVICE HELP LINE

For technical assistance, call Moldex Technical Service Department at +1 (310) 837-6500 ext. 512/550 or +1 (800) 421-0668 ext. 512/550.


MOLDEX-METRIC, INC.

10111 W. Jefferson Blvd. Culver City, CA 90232
 TEL: +1 (800) 421-0668 or +1 (310) 837-6500 FAX +1 (310) 837-9563
 www.moldex.com e-mail: sales@moldex.com


CANADIAN CUSTOMER SERVICE

TEL: +1 (310) 837-6500, Ext. 517
 FAX: +1 (310) 837-9563

Moldex, Ideas that wear well, Softspun, Dura-Mesh, Flexwing, Ventex and EZ-ON are registered trademarks of Moldex-Metric, Inc. The mesh design is a trademark of Moldex-Metric, Inc. U.S. Patents #D285,374, #4,850,347, #6,571,797, #D389,239, and #6,047,698. Made in the U.S.A. 9702-712 REV E 06/09



Moldex-Metric, Inc.
 CULVER CITY, CALIFORNIA, 90232
 +1 (800) 421-0668



THESE RESPIRATORS ARE APPROVED ONLY IN THE FOLLOWING CONFIGURATIONS:

TC	PROTECTION	RESPIRATOR					CAUTIONS AND LIMITATIONS*
		NOSE	EYE	FACE	HEAD	HAIR	
844-3743	N95	X	X				A B C J M N O P
844-3846	N95			X	X		A B C J M N O P

1. Protection
 N95-Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply.

2. Cautions and Limitations
 A- Not for use in atmospheres containing less than 19.5 percent oxygen.
 B- Not for use in atmospheres immediately dangerous to life or health.
 C- Do not exceed maximum use concentrations established by regulatory standards.
 J- Failure to properly use and maintain this product could result in injury or death.
 M- All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
 N- Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
 O- Refer to User's Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
 P- NIOSH does not evaluate respirators as surgical masks.

PACKAGING

EZ 22, Medium/Large
 10 Respirators per Bag
 20 Bags per Case

EZ 23, Medium/Large
 10 Respirators per Bag
 10 Lockers per Case

EZ 22S, Small
 10 Respirators per Bag
 20 Bags per Case

EZ 23S, Small
 10 Respirators per Bag
 10 Respirators per Case

LIMITED WARRANTY IMPORTANT NOTICE TO PURCHASER

This limited warranty is made in lieu of the warranties of merchantability, fitness for particular purposes and all other warranties, express or implied. There are no other warranties which extend beyond the description on the face hereof. The physical standards and specifications of Moldex will be met by products sold. **Exclusive Remedies:** damages for the breach of this limited warranty are limited to the replacement of such quantity of Moldex products proved to be defectively manufactured. Except as provided above, Moldex shall not be liable or responsible for any loss, damages, or liability, direct, indirect, incidental, special or consequential, arising out of a sale, use or misuse, or the inability to use products by the user.

WARNING TO USER

1. Follow all instructions and warnings on the use of this respirator and wear during all times of exposure. Failure to do so will reduce respirator effectiveness, wearer protection and may result in sickness or death.
2. Before use, the user must first be trained by the employer in proper respirator use, in accordance with applicable safety and health standards, for the contaminant and exposure level in the assigned work area.
3. The particulates and other contaminants which can be dangerous to your health include those which you cannot see, taste, or smell.
4. This product has not been sold with warnings or use instructions for personnel involved in healthcare or related situations, where there may be the possibility of contact with disease or biological hazards. If you are considering such uses, first call the Moldex Technical Service Dept., +1 (310) 837-6500 ext. 512/550 or +1 (800) 421-0668 ext. 512/550.

DISTRIBUTED BY:

