

# Disposable Respiratory Protection for Fire and EMS Services



**1500 N95**



**The EMS Necessity**



**2730N100**



**Overhaul Respiratory Protection**



**2360P100**



**Global Pandemic Preparedness**



**Soft foam nose cushion eliminates pressure points**

**Easy donning - no adjustment of noseband**

**Surgical N95 NIOSH Certified Respirator**  
(Surgical N95 applies to 1500 N95 Series only)



**Softspun® inner lining for added comfort**

**Dura-Mesh® shell resists collapsing**



**1500 N95 Series for EMS Transport**

- Surgical N95 NIOSH Certified Respirator.
- Meets CDC guidelines for protection against TB.
- Dura-Mesh® shell resists collapsing in regular use and even under hot and humid conditions.
- Contour molded nosebridge for extra comfort and easier donning.
- Latex Free.
- Highest ASTM level of fluid resistance on the market (160mm Hg).
- Available in 5 color coded sizes.



**2730N100 Series High Efficiency Filtration for Overhaul Environments**

- Dura-Mesh® shell resists collapsing in regular use and even under hot and humid conditions.
- NIOSH certified to have a filter efficiency of 99.97% effective against non-oil based particulates.
- Exclusive Ventex® valve lets the hot air out faster so wearers breathe easier and stay cooler longer.
- Available in two sizes and two types of packaging.



**2360P100 Pandemic Preparedness and High Efficiency Filtration**

- Dura-Mesh® shell resists collapsing in regular use and even under hot and humid conditions.
- NIOSH certified to have a filter efficiency of 99.97% effective against oil and non-oil based particulates.
- Adjustable cloth straps allow a comfortable fit with a simple tug on the straps.

<u>Part No.</u>	<u>Size</u>	<u>Quantity</u>
1510	Extra Small	20 Masks per Box
1511	Small	8 Boxes per Case
1512	Medium	
1513	Large	
1517	Low Profile	
2730N100	Medium/Large	5 Masks per Box
2731N100	Small	6 Boxes per Case
2730AN100	Medium/Large	Individually Packaged
2731AN100	Small	20 Masks per Case
2360P100	Medium/Large	5 Masks per Bag 4 Bags per Case

**Distributed By:**



**Moldex-Metric, Inc. • 10111 West Jefferson Blvd. • Culver City, CA 90232**  
Phone: (800) 421-0668 • (310) 837-6500 Ext. 512/550 • Fax: (310) 837-9563  
Email: sales@moldex.com • www.moldex.com

